

Breast Mass Analyzer (BMA)

Background

Breast cancer represents the most common oncological disease in women worldwide, and a multimodality imaging approach is becoming increasingly important to detect and monitor early-stage cancer and to initiate timely and appropriate treatment in order to reduce mortality. According to WHO data, in 2020, there were 2.3 million women diagnosed with breast cancer and 685,000 deaths globally^[1].

Technological developments over recent years have made breast ultrasound (US) an essential component for the breast cancer care, together with mammography/tomosynthesis and magnetic resonance imaging (MRI). Mammography is the gold standard screening technique for breast cancer^[2]. Ultrasound is not typically used at this stage but appears to be especially helpful in mass characterization where mammogram performance is limited and can be difficult to see abnormal features, as well as when a better look at a suspicious area, previously identified on mammography, is required^[3].

Other advantages of US include its wide availability, good tolerance by patients, lack of ionizing radiation and relative low cost^[4].

Moreover, US provides an optimal, cheap, and comfortable guidance for performing percutaneous needle biopsy, and, in particular, large-core needle biopsy (US-LCNB), the gold standard in diagnosis in case of suspect of breast cancer. This safe and cost-effective procedure has almost entirely replaced surgical biopsy, thus limiting costs, patients discomfort and possible complications^[5]. However, it remains a semi-invasive procedure associated with an anticipatory emotional distress which could be avoided, where not necessary.

BI-RADS® (Breast Imaging-Reporting and Data System) provides a widely accepted reporting schema for imaging of the breast and characterization of findings. This system is designed to standardize breast imaging reporting and to reduce confusion in breast imaging interpretation, thus facilitating outcome monitoring and quality assessment. Results are sorted into categories numbered 0 through 6, by increasing the order of likelihood of findings ^[6,7].

Category		Management	Likelihood of cancer
0	Need additional imaging or prior examinations	Recall for additional imaging and/or await prior examinations	N/A
1	Negative	Routine screening	Essentially 0%
2	Benign	Routine screening	Essentially 0%
3	Probably Benign	Short interval-follow-up (6 months) or continued	>0% but ≤2%
4	Suspicious	Tissue sampling	4a. low suspicion for malignancy (>2% to ≤10%); 4b. moderate suspicion for malignancy (>10% to ≤50%); 4c. high suspicion for malignancy (>50% to <95%)
5	Highly suggestive of malignancy	Tissue sampling	≥95%
6	Known biopsy-proven	Surgical excision when clinical appropriate	N/A

Table 1: BI-RADS® final assessment category.

BI-RADS® category 3 represents a finding with a very low ($>0\%$ but $\leq 2\%$) probability of being cancer, that is not expected to change over time. A follow-up with repeat imaging in 3 to 6 months and regularly after that until the finding is known to be stable (usually at least 2 years) is usually prescribed. Studies aiming at reducing the US-guided biopsy rate should take particular attention to this category.

US images contain a high number of quantitative features, known as radiomics, which can be associated with a higher or lower probability of malignancy and may thus help to better characterize breast lesions^[8-10]. These features include not only shape, margins, and spatial orientation of lesions, but also their signal intensity and heterogeneity, the quantification of which is difficult for the human reader. Recent technological advances in the field of artificial

intelligence applied to image analysis allow the development of machine learning models predicting different class of risks based on the wide spectrum of information found in a US image^[11-12].

TRACE4BUS is a machine learning model based on radiomics features belonging to different families, such as morphology, intensity-based statistics, intensity histogram, grey-level co-occurrence matrix, grey-level run length matrix, grey-level size zone matrix^[13-14], neighborhood grey tone difference matrix, grey-level distance zone matrix, and neighboring grey-level dependence matrix. Their definition, computation, and nomenclature are compliant with the International Biomarker Standardization Initiative (IBSI) guidelines^[15].

Rank	Feature Family	Feature Name
1	Morphology	Perimeter-to-area ratio **
2	Morphology	Maximum diameter **
3	Morphology	Compactness **
4	Morphology	Acircularity **
5	Morphology	Perimeter **
6	Morphology	Area **
7	Morphology	Center of mass shift **
8	Morphology	Circularity *
9	Neighborhood grey tone difference matrix	Strength **
10	Neighborhood grey tone difference matrix	Coarseness **
11	Neighborhood grey tone difference matrix	Contrast
12	Neighborhood grey tone difference matrix	Busyness *
13	Grey-level size zone matrix	Zone size non-uniformity **
14	Grey-level size zone matrix	Grey-level non-uniformity glszm **
15	Neighboring grey-level dependence matrix	Dependence count non-uniformity **
16	Neighboring grey-level dependence matrix	Low-dependence low-grey-level emphasis
17	Grey-level run length matrix	Grey-level non-uniformity
18	Grey-level run length matrix	Run length non-uniformity
19	Intensity-based statistics	Minimum
20	Intensity-based statistics	Energy
21	Intensity-based statistics	Variance
22	Intensity-based statistics	Quartile coefficient
23	Intensity-based statistics	10th percentile
24	Intensity histogram	10th percentile
25	Grey-level co-occurrence matrix	First measure of information correlation

* Denotes statistical significance at 0.05 (adjusted with Bonferroni-Holm correction).

** denotes a statistical significance at 0.005 (adjusted with Bonferroni-Holm correction).

Table 2: Top 25 most relevant predictors sorted in descending order of relevance.

TRACE4BUS: clinical validation

The development and the clinical validation of TRACE4BUS is described in detail in the paper by Interlenghi M. et al. "A Machine Learning Ensemble Based on Radiomics to Predict BI-RADS® Category and Reduce the Biopsy Rate of Ultrasound-Detected Suspicious Breast Masses"^[12]. Esaote systems were part of the pool of ultrasound devices that were used for the development and validation of the algorithm.



Article

A Machine Learning Ensemble Based on Radiomics to Predict BI-RADS Category and Reduce the Biopsy Rate of Ultrasound-Detected Suspicious Breast Masses

Matteo Interlenghi^{1,†}, Christian Salvatore^{1,2,†}, Veronica Magni³, Gabriele Caldara², Elia Schiavon¹, Andrea Cozzi³, Simone Schiaffino⁴, Luca Alessandro Carbonaro^{5,6}, Isabella Castiglioni^{7,8,*} and Francesco Sardanelli^{3,4}

- ¹ DeepTrace Technologies S.R.L., Via Conservatorio 17, 20122 Milano, Italy; interlenghi@deeptech.com (M.I.); salvatore@deeptech.com (C.S.); schiavon@deeptech.com (E.S.)
 - ² Department of Science, Technology and Society, Scuola Universitaria IUSS, Istituto Universitario di Studi Superiori, Piazza della Vittoria 15, 27100 Pavia, Italy; gabriele.caldara@iusspavia.it
 - ³ Department of Biomedical Sciences for Health, Università degli Studi di Milano, Via Luigi Mangiagalli 31, 20133 Milano, Italy; veronica.magni@unimi.it (V.M.); andrea.cozzi1@unimi.it (A.C.); francesco.sardanelli@unimi.it (F.S.)
 - ⁴ Unit of Radiology, IRCCS Policlinico San Donato, Via Rodolfo Morandi 30, 20097 San Donato Milanese, Italy; simone.schiaffino@grupposandonato.it
 - ⁵ Department of Radiology, ASST Grande Ospedale Metropolitano Niguarda, Piazza dell'Ospedale Maggiore 3, 20162 Milano, Italy; luca.carbonaro@unimi.it
 - ⁶ Department of Oncology and Hemato-Oncology, Università degli Studi di Milano, Via Festa del Perdono 7, 20122 Milan, Italy
 - ⁷ Institute of Biomedical Imaging and Physiology, Consiglio Nazionale delle Ricerche, Via Fratelli Cervi 93, 20090 Segrate, Italy
 - ⁸ Department of Physics, Università degli Studi di Milano-Bicocca, Piazza Della Scienza 3, 20126 Milano, Italy
- * Correspondence: isabella.castiglioni@unimib.it
- † These authors contributed equally to this work.

Abstract: We developed a machine learning model based on radiomics to predict the BI-RADS category of ultrasound-detected suspicious breast lesions and support medical decision-making towards short-interval follow-up versus tissue sampling. From a retrospective 2015–2019 series of ultrasound-guided core needle biopsies performed by four board-certified breast radiologists using six ultrasound systems from three vendors, we collected 821 images of 834 suspicious breast masses from 819 patients, 404 malignant and 430 benign according to histopathology. A balanced image set of biopsy-proven benign ($n = 299$) and malignant ($n = 299$) lesions was used for training and cross-validation of ensembles of machine learning algorithms supervised during learning by histopathological diagnosis as a reference standard. Based on a majority vote (over 80% of the votes to have a valid prediction of benign lesion), an ensemble of support vector machines showed an ability to reduce the biopsy rate of benign lesions by 15% to 18%, always keeping a sensitivity over 94%, when externally tested on 236 images from two image sets: (1) 123 lesions (51 malignant and 72 benign) obtained from two ultrasound systems used for training and from a different one, resulting in a positive predictive value (PPV) of 45.9% (95% confidence interval 36.3–55.7%) versus a radiologists' PPV of 41.5% ($p < 0.005$), combined with a 98.0% sensitivity (89.6–99.9%); (2) 113 lesions (54 malignant and 59 benign) obtained from two ultrasound systems from vendors different from those used for training, resulting into a 50.5% PPV (40.4–60.6%) versus a radiologists' PPV of 47.8% ($p < 0.005$), combined with a 94.4% sensitivity (84.6–98.8%). Errors in BI-RADS 3 category (i.e., assigned by the model as BI-RADS 4) were 0.8% and 2.7% in the *Testing set I* and *II*, respectively. The board-certified breast radiologist accepted the BI-RADS classes assigned by the model in 114 masses (92.7%) and modified the BI-RADS classes of 9 breast masses (7.3%). In six of nine cases, the model performed better than the radiologist did, since it assigned a BI-RADS 3 classification to histopathology-confirmed benign masses that were classified as BI-RADS 4 by the radiologist.



Citation: Interlenghi, M.; Salvatore, C.; Magni, V.; Caldara, G.; Schiavon, E.; Cozzi, A.; Schiaffino, S.; Carbonaro, L.A.; Castiglioni, I.; Sardanelli, F. A Machine Learning Ensemble Based on Radiomics to Predict BI-RADS Category and Reduce the Biopsy Rate of Ultrasound-Detected Suspicious Breast Masses. *Diagnostics* **2022**, *12*, 187. <https://doi.org/10.3390/diagnostics12010187>

Academic Editor: Dechang Chen

Received: 2 December 2021

Accepted: 30 January 2022

Published: 13 January 2022

Publisher's Note: MDPI stays neutral with regard to jurisdictional claims in published maps and institutional affiliations.



Copyright: © 2022 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (<https://creativecommons.org/licenses/by/4.0/>).

Diagnostics **2022**, *12*, 187. <https://doi.org/10.3390/diagnostics12010187>

<https://www.mdpi.com/journal/diagnostics>

Fig. 1: Clinical study front page.

Augmented Insight™, an Esaote smart solution applied to breast ultrasound imaging

Augmented Insight™ embeds all the solutions powered by Artificial Intelligence, developed by the Research & Development department at Esaote according to a multidisciplinary approach, across different modalities such as Ultrasound, MRI and Ebit. Based on data-driven machine learning, Augmented Insight™ is designed to simplify the workflow in repetitive gestures, complex procedures, measurements, or lesions analysis. Esaote, particularly active in the breast imaging field and in constant search of innovative solutions, has recently developed many Augmented Insight™ solutions specifically for breast US imaging, including BreastNav™MRI, eDetect and Breast Mass Analyzer technologies.

Esaote Breast Mass Analyzer Technology

Esaote Breast Mass Analyzer (BMA) has embedded the TRACE4BUS software by Deep Trace Technologies. BMA Augmented Insight™ technology is a radiomics-based machine learning model to predict the BI-RADS® category of ultrasound-detected suspicious breast lesions and support medical decision-making towards BI-RADS® 3 versus BI-RADS® 4–5.

Therefore, this technology is indicated for all physicians with experience in breast ultrasound seeking support for reporting ultrasound-detected suspicious breast masses, assigning the diagnostic category to such masses according to BI-RADS® classification, and, consequently, making an informed decision whether or not to send such masses for ultrasound-guided needle biopsy.

During the ultrasound examination, as soon as the user identifies a suspected lesion, freezes the image and manually contours the mass using the trackball.

The software, supported by A.I., automatically processes the mass, and within a few seconds provides its characteristics (i.e. shape and orientation) and BI-RADS® category (BI-RADS® category calculated by BMA is on 3 levels: 3, 4 and 5).

The user can subsequently accept the proposed value, modify, or refine it on categories 4a, 4b and 4c, directly from the touchscreen. If approved, all the information are automatically included in the worksheet environment and in the report.

Innovation and research

This innovative BMA tool further expands the scope of the company's breast ultrasound technologies dedicated to the breast. It has proved its viability as a second opinion on BI-RADS® classification as well as its potential clinical value in reducing the biopsy rate.

The underlying philosophy of research and development in Esaote is based on open innovation involving external expertise in a valuable co-operative network that includes research centers and highly experienced physicians and their patients.




Fig. 2: BMA procedure: step 1 - mass processing; step 2 - BI-RADS® categorization.

The report, a dedicated template (that can be edited directly from the touchscreen) produced according to the BI-RADS® scoring and categorization, provides information on the localization of the lesion, with a dedicated body mark to visually locate the lesion in the correct quadrant, as well as details on the mass itself (shape, orientation, and echo pattern) and its margins (circumscribed, lobulated and spiculated).

To enhance the quality of lesion characterization, additional technologies such as microV, the Esaote tool for microvascularization visualization, and QElaXto 2D, the Esaote 2D shear wave technique for stiffness evaluation of the breast tissue, are available on the system.



Breast					
R Mass			L Mass		
R Mass 1			L Mass 1		
Length	6.8	mm	Length	3.3	mm
Height	8.6	mm	Height	6.7	mm
Width	10.7	mm	Width	1.9	mm
Volume	0.328	ml	Volume	0.022	ml
Skin Depth	5.8	mm	Skin Depth	20.7	mm
Nipple Dist	10.0	mm	Nipple Dist	23.0	mm
Location			Location		
O'clock	---		O'clock	---	
Region	Nipple		Region	Areolar	
Quadrants	Upper Inner		Quadrants	Upper Outer	
Profile	Posterior		Profile	Middle	
Probe Angle	60		Probe Angle	120	
Masses			Masses		
Shape	Oval		Shape	Oval	
Orientation	Parallel		Orientation	Parallel	
Echo Pattern	Anechoic		Echo Pattern	Hypoechoic	
Posterior Features	No posterior features		Posterior Features	No posterior features	
Margin			Margin		
Circumscribed	Yes		Circumscribed	Yes	
Not circumscribed - Indistinct	No		Not circumscribed - Indistinct	No	
Not circumscribed - Angular	No		Not circumscribed - Angular	No	
Not circumscribed - Microlobulated	No		Not circumscribed - Microlobulated	No	
Not circumscribed - Spiculated	No		Not circumscribed - Spiculated	No	
Elasticity Assess	---		Elasticity Assess	---	
BI-RADS Category	3 - Probably Benign		BI-RADS Category	3 - Probably Benign	
					

Fig. 3: BI-RADS® report.

References

1. Mangone L, Bisceglia I, Michiara M, Musolino A, Mazzoleni G, Caldarella A, Minerba S, Cascone G, Bella F, Dinaro Y, Pau L, Pinto C. Breast Cancer in Italy: Stage and Region Distribution. *Breast Cancer* (Dove Med Press). 2022 Apr 29;14:125-131. doi: 10.2147/BCTT.S360244. PMID: 35515355; PMCID: PMC9064450;
2. Miller RG. Breast cancer screening: can we talk? *J Gen Intern Med*. 2001 Mar;16(3):206-7. doi: 10.1111/j.1525-1497.2001.10119.x. PMID: 11318918; PMCID: PMC1495188;
3. American College of Radiology. ACR BI-RADS ATLAS – Mammography. Reporting System, 2013. Accessed at <https://www.acr.org/-/media/ACR/Files/RADS/BI-RADS/Mammography-Reporting.pdf> on September 30, 2021.
4. Jaime Geisel, Madhavi Raghu, Regina Hooley, The Role of Ultrasound in Breast Cancer Screening: The Case for and Against Ultrasound, *Seminars in Ultrasound, CT and MRI*, Volume 39, Issue 1, 2018, Pages 25-34, ISSN 0887-2171, <https://doi.org/10.1053/j.sult.2017.09.006>;
5. Bick, U.; Trimboli, R.M.; Athanasiou, A.; Balleyguier, C.; Baltzer, P.A.T.; Bernathova, M.; Borbély, K.; Brkljacic, B.; Carbonaro, L.A.; Clauser, P.; et al. Image-guided breast biopsy and localisation: Recommendations for information to women and referring physicians by the European Society of Breast Imaging. *Insights Imaging* 2020, 11, 12.
6. Breast Imaging Reporting and Data System? (BI-RADS?) Atlas.
7. Breast Imaging Reporting and Data System, Inter- and Intraobserver Variability in Feature Analysis and Final Assessment by Wendie A. Berg et al. Department of Radiology, University of Maryland School of Medicine, 22 S. Greene St., Baltimore, *AJR* 2000; 174:1769-1777.
8. Tagliafico, A.S.; Piana, M.; Schenone, D.; Lai, R.; Massone, A.M.; Houssami, N. Overview of radiomics in breast cancer diagnosis and prognostication. *Breast* 2020, 49, 74-80;
9. Stavros, A.T.; Thickman, D.; Rapp, C.L.; Dennis, M.A.; Parker, S.H.; Sisney, G.A. Solid breast nodules: Use of sonography to distinguish between benign and malignant lesions. *Radiology* 1995, 196, 123-134;
10. Bradley J. Erickson, Panagiotis Korfiatis, Zeynettin Akkus, Timothy L. Kline, *RadioGraphics*, Feb 17 2017 <https://doi.org/10.1148/rg.2017160130>;
11. Bitencourt, A.; Daimiel Naranjo, I.; Lo Gullo, R.; Rossi Saccarelli, C.; Pinker, K. Alenhanced breast imaging: Where are we and where are we heading? *Eur. J. Radiol.* 2021, 142, 109882. [Google Scholar] [CrossRef];
12. Hu, Q.; Giger, M.L. Clinical Artificial Intelligence Applications: Breast Imaging. *Radiol. Clin. N. Am.* 2021, 59, 1027-1043. <https://doi.org/10.1016/j.breast.2019.10.018>;
13. TRACE4. Available online: http://www.deeptracetech.com/temp/TechnicalSheet__TRACE4.pdf (accessed on 7 January 2022);
14. Interlenghi M, Salvatore C, Magni V, Caldara G, Schiavon E, Cozzi A, Schiaffino S, Carbonaro LA, Castiglioni I, Sardanelli F. A Machine Learning Ensemble Based on Radiomics to Predict BI-RADS® Category and Reduce the Biopsy Rate of Ultrasound-Detected Suspicious Breast Masses. *Diagnostics* (Basel). 2022 Jan 13;12(1):187. doi: 10.3390/diagnostics12010187. PMID: 35054354; PMCID: PMC8774734;
15. Zwanenburg, A.; Vallières, M.; Abdalah, M.A.; Aerts, H.J.W.L.; Andrearczyk, V.; Apte, A.; Ashrafinia, S.; Bakas, S.; Beukinga, R.J.; Boellaard, R.; et al. The Image Biomarker Standardization Initiative: Standardized Quantitative Radiomics for High Throughput Image-based Phenotyping. *Radiology* 2020, 295, 328-338

Notes



160000514 ver 01



Esaote S.p.A. - sole-shareholder company

Via Enrico Melen 77, 16152 Genova, ITALY, Tel. +39 010 6547 1, Fax +39 010 6547 275,
info@esaote.com - www.esaote.com

BMA is powered by Deep Trace Technology. BI-RADS® incorporates the Breast Imaging Reporting and Data System ATLAS of the American College of Radiology, Copyright 1992, 1993, 1995, 1998, 2003, and 2013. The developer of this product is independently owned and operated, and is not an affiliate of the American College of Radiology. The American College of Radiology is not responsible for the contents or operation of this product or its associated software, and expressly disclaims any and all warranties and liabilities, expressed or implied, in connection therewith. Technology and features are system/configuration dependent. Specifications subject to change without notice. Information might refer to products or modalities not yet approved in all countries. Product images are for illustrative purposes only. For further details, please contact your Esaote sales representative.

Please visit us online
for more information



Italian design