


iQProbes

Women's Health Ultrasound:
Advanced solutions
in probe technology



New iQProbe technology for high-quality imaging in Women's Health

Ramona De Luca, PhD, Esaote Ultrasound Probes R&D (Florence, Italy)

Martina Cereseto, Esaote Global Marketing (Genoa, Italy)



The C 2-9 is a high frequency convex array with a specific configuration responsible for perfect characterization of foetus tiny details in the first trimester ultrasound scan. It also provides anatomical scan in the second trimester with great clarity though the whole image field of view.



The E 3-12 is an endovaginal endfire transducer with enhanced sensitivity and high resolution. It offers high clinical confidence in assessing anomalies and pathologies of the uterus and the ovaries and aids the early pregnancy imaging.



Esaote offers a versatile portfolio of probes for high-quality imaging of Women's Healthcare

Introduction

As a leader in design and manufacturing of ultrasound probes, Esaote is committed to R&D investment and innovations, and gathers feedback from users to continually improve ergonomics, reliability, and image quality. Esaote's ultrasound systems and probes are designed to make easy and confident diagnoses across a variety of clinical applications, without compromising the comfort of clinician or patient. Recent introduction of new ultrasound transducer architectures with advanced materials and innovative designs is providing high-quality imaging for Women's Healthcare. The new C 2-9 and E 3-12 probes incorporate several innovative technologies that translate into exceptional resolution, penetration, and overall uniformity in all imaging modes.

Probe description

Ultrasound probes (Figure 1) are critical to achieve exceptional imaging performance and repeatability. Their design requires development of advanced materials in order to improve transducer operation efficiency and deliver remarkable image quality (IQ) performances. Recent technology embedded in the C 2-9 and E 3-12 applies a combination of advanced materials together with specifically designed transducer geometry. An acoustic lens material is used that minimizes reverberations and enhances image contrast resolution. Esaote's innovative backing block increases ultrasound energy transmitted into the patient body, while maintaining very wide bandwidth (Figure 2). This directly translates into enhanced image sensitivity, higher resolution, useful penetration, and overall clarity from near field to far field. Automation and aggressive reprocessing (cleaning, disinfection and sterilization) trends require improved probe reliability. In response, Esaote developed a special protective layer placed beneath the acoustic lens which effectively protects the transducer and prevents fluid ingress. With these innovations, Esaote is delivering a new level of diagnostic confidence and accuracy in all major modes – whether fundamental imaging, Doppler, or tissue harmonic imaging – and across the entire range of clinical applications.

Figure 1 - Medical ultrasound transducers consist of array of small piezoelectric elements attached to three other structures: matching layers, backing block and acoustic lens.





Figure 2 – New generation probes provide new levels of IQ performances in different types of modes – fundamental imaging, Doppler, and harmonic imaging.



Clinical use

1st trimester pregnancy



1st trimester fetus intracranial structures



2nd trimester fetus diaphragm



2nd trimester fetal heart



3rd trimester fetal heart B-mode and CFM



2nd trimester fetus kidney with CFM



1st trimester twins




2nd trimester circle of Willis with microV




2nd trimester fetus renal arteries with microV




Ovary perfusion with microV




ElaXto on uterus



Ovary with Power Doppler




Early pregnancy endocavity approach



Early pregnancy fetus and gestational sac




Early pregnancy dual mode



Ovary and follicles



Healthy uterus



Ovary B-mode



Esaote S.p.A. - sole-shareholder company

Via Enrico Meloni 77, 16152 Genova, ITALY, Tel. +39 010 6547 1, Fax +39 010 6547 275, info@esaote.com

Technology and features are system/configuration dependent. Specifications subject to change without notice. Information might refer to products or modalities not yet approved in all countries. Product images are for illustrative purposes only. For further details, please contact your Esaote sales representative.