



<b>IHE Integration Statement</b>		<b>Date:</b> Mar. 5 <sup>th</sup> , 2024
<b>Vendor</b>	<b>Product Name</b>	<b>SW Build</b>
ESAOTE S.p.A.	MyLab™ X90, X9, 9, X8, X75, A70, A50, X7, X6, X5, Omega and Sigma Platforms, XPro80, XPro30 *	F12XXXX
These products are intended to implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>
Scheduled Workflow	Acquisition Modality	Broad Worklist Query
Patient Information Reconciliation	Acquisition Modality	none
Portable Data for Imaging	Portable Media Creator	none
Evidence Documents	Acquisition Modality	none
Echocardiography Workflow	Acquisition Modality	Broad Worklist Query Stress Echo
Cardiology Evidence Documents	Acquisition Modality	Echo Evidence
Consistent Time	Time Client	none
Internet address for vendor's IHE information: <a href="http://www.esaote.com">http://www.esaote.com</a>		
Vendor's e-mail address for information about IHE: <a href="mailto:ihe@esaote.com">ihe@esaote.com</a>		
<b>Links to Standards Conformance Statements for the Implementation</b>		
<b>DICOM</b> ®	<a href="http://www.esaote.com/dicom.htm">http://www.esaote.com/dicom.htm</a>	
<b>Links to general information on IHE</b>		
In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">www.ihe-j.org</a>

\* all with DICOM® option.

ESAOTE S.p.A.  
via Melen 77  
16152 Genova (Italy)

© Copyright ESAOTE S.p.A. 2024. All rights reserved.

\* DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.